

DUST assay in human cells

YS Yilun Sun JN John L. Nitiss YP Yves Pommier

Updated date: Oct 8, 2021

 An abbreviated version of this protocol was published in Science Advances in Nov 2020

A conserved SUMO pathway repairs topoisomerase DNA-protein cross-links by engaging ubiquitin-mediated proteasomal degradation

DOI: 10.1126/sciadv.aba6290

Detailed protocol

1. Seed 5×10^5 HEK293 cells in 35-mm dish per sample the day before treatment.
2. Next day, add topoisomerase inhibitors of desired concentrations to the cells.
3. Collect cells after treatment by washing with $1 \times$ PBS then lysing with 600 μ l of DNAzol (Invitrogen) containing $1 \times$ protease cocktail inhibitor and 1 mM DTT for 10 min on ice.
4. Add 300 μ l 200 proof ethanol to the dishes and place the dishes on a shaker at 4 °C till white aggregates (dehydrated DNA) become visible.
5. Transfer the lysates to 1.5 mL eppendorf tube for centrifuge at 15,000 rpm at 4 °C.
6. Remove the supernatant and wash the nucleic acid pellet at the bottom of the tube with 75% ethanol. Remove the ethanol and resuspend the pellet in 200 μ l of TE buffer,
7. Heat the nucleic acid sample at 37°C for 15 min, followed by shearing with sonication (40% output for 10-s pulse and 10-s rest for four times).
8. Centrifuge the samples at 15,000 rpm for 5 min at 4°C, collect the supernatant and treat with RNase A (100 μ g/ml) for 1 hour at 37°C to remove RNA contamination (**the RNA removal step is optional if the yield of nucleic acid is low**), followed by addition of 1:10 volume of 3 M sodium acetate sodium acetate and 2.5 volume of 200 proof ethanol.
9. After 20 min 15,000 rpm centrifugation, recover and resuspend the DNA pellet in 100 μ l of ddH₂O.
10. Take 1 μ l per sample for spectrophotometric measurement of absorbance at 260 nm to quantitate DNA content (NanoDrop).
11. DNA (10 μ g) from each sample was digested with 50 units of micrococcal nuclease (100 units/ μ l; Thermo Fisher Scientific) in the presence of 5 mM CaCl₂ at 37°C for 30 min.
12. Add 4 \times SDS sample loading buffer to the digested samples, followed by gel electrophoresis on 4 to 15% precast polyacrylamide gel (Bio-Rad) for immunodetection for total TOP-DPCs (Anti-TOP1, mouse monoclonal BD Biosciences Cat# 556597, 1:500; Anti-TOP2 α , mouse monoclonal Millipore Cat# MAB4197, 1:500; Anti-TOP2 β , mouse monoclonal BD Biosciences Cat# 611492, 1:500), SUMO-TOP-DPCs (Anti-SUMO-1, rabbit monoclonal Cell Signaling Cat# 4940, 1:250; Anti-SUMO-2/3, rabbit monoclonal Cell Signaling Cat# 4971, 1:250, and Ub-TOP1-DPCs (Anti-ubiquitin, mouse monoclonal Santa Cruz Cat# sc-8017, 1:100;)
13. 2 μ g of each sample (without micrococcal nuclease digestion) was subjected to slot-blot for immunoblotting with anti-double-stranded DNA antibody (Anti-dsDNA, mouse monoclonal Abcam Cat# ab27156, 1:5000) as a loading control to verify that amounts of DNA were digested with micrococcal nuclease.

How to cite: (Readers should cite both the Bio-protocol preprint and the original research article where this protocol was used)

1. Sun, Y. , Nitiss, J. and Pommier, Y. (2021). DUST assay in human cells. Bio-protocol Preprint. bio-protocol.org/prep1399.
2. Sun, Y., Jenkins, L. M. M., Su, Y. P., Nitiss, K. C., Nitiss, J. L. and Pommier, Y. (2020). A conserved SUMO pathway repairs topoisomerase DNA-protein cross-links by engaging ubiquitin-mediated proteasomal degradation . Science Advances 6(46). DOI: [10.1126/sciadv.aba6290](https://doi.org/10.1126/sciadv.aba6290)

Copyright: Content may be subjected to copyright.